

Delta-Lite[®] Plus

The Trusted Rigid Solution



Unique Features

- Very smooth surface
- Excellent rigid support
- High lamination
- New resin formulation - use exam glove of your choice
- Weight bearing in approximately 20 minutes
- 3-5 minute set time

Delta-Lite® Plus

The Trusted Rigid Solution

Your new standard of casting.

Delta-Lite® Plus is our new fiberglass cast tape that combines a fiberglass substrate and advanced water activated polyurethane resin system which provides tack free properties. The tack free properties allow ease of application and molding using standard examination gloves.



Features and Benefits

Very Smooth surface - Less abrasive finished cast enhances patient comfort while it reduces clothing snags

Excellent rigid support (immobilisation) - Creates long-lasting, light weight cast

High lamination - Creates a strong bond between layers ensuring a strong finished cast

New resin formulation - Use exam glove of your choice. Improves unrolling tension improving ease of cast application

Good conformability - Makes casting application easier with change of anatomy

Outstanding moulding - Produces an anatomically correct cast

Ordering Information

Order Number	Size	Color	Unit of Measure
73458-00	1 in. x 2 yds.	WHITE	10 rolls / box
73458-01	2 in. x 4 yds.		10 rolls / box
73458-02	3 in. x 4 yds.		10 rolls / box
73458-03	4 in. x 4 yds.		10 rolls / box
73458-04	5 in. x 4 yds.		10 rolls / box
73458-45	2 in. x 4 yds.	BLACK	10 rolls / box
73458-46	3 in. x 4 yds.		10 rolls / box
73458-47	4 in. x 4 yds.		10 rolls / box
73458-35	2 in. x 4 yds.	LIGHT BLUE	10 rolls / box
73458-36	3 in. x 4 yds.		10 rolls / box
73458-37	4 in. x 4 yds.		10 rolls / box
73458-38	5 in. x 4 yds.		10 rolls / box
73458-60	2 in. x 4 yds.	PURPLE	10 rolls / box
73458-61	3 in. x 4 yds.		10 rolls / box
73458-62	4 in. x 4 yds.		10 rolls / box
73458-55	2 in. x 4 yds.	PINK	10 rolls / box
73458-56	3 in. x 4 yds.		10 rolls / box
73458-57	4 in. x 4 yds.		10 rolls / box
73458-30	2 in. x 4 yds.	DEEP RED	10 rolls / box
73458-31	3 in. x 4 yds.		10 rolls / box
73458-32	4 in. x 4 yds.		10 rolls / box
73458-20	2 in. x 4 yds.	DEEP BLUE	10 rolls / box
73458-21	3 in. x 4 yds.		10 rolls / box
73458-22	4 in. x 4 yds.		10 rolls / box
73458-23	5 in. x 4 yds.		10 rolls / box
73458-25	2 in. x 4 yds.	DEEP GREEN	10 rolls / box
73458-26	3 in. x 4 yds.		10 rolls / box
73458-27	4 in. x 4 yds.		10 rolls / box
73458-28	5 in. x 4 yds.		10 rolls / box
73458-65	2 in. x 4 yds.	YELLOW	10 rolls / box
73458-66	3 in. x 4 yds.		10 rolls / box
73458-67	4 in. x 4 yds.		10 rolls / box
73458-50	2 in. x 4 yds.	MIXED PACK*	10 rolls / box
73458-51	3 in. x 4 yds.		10 rolls / box
73458-52	4 in. x 4 yds.		10 rolls / box

*Mixed pack: Contains 2 Deep Blue, 2 Light Blue, 2 Deep Red, 2 Purple and 2 Pink.

Colors and availability subject to change. Please refer to the product label and / or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
 Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
 To order toll-free: BSN medical (+1) 800 552 1157