


# Leukomed® sterile wound dressings




**Ultra-thin** for high vapour transmission





The **pattern** for breathable film dressings



**Ease** of application

| Leukomed®  | Order No. | HCPCS | Dressing Size  | Wound Pad Size    | Dressings/Box | Shipping Unit |
|--|-----------|-------|----------------|-------------------|---------------|---------------|
|  | 7238001   | A6254 | 3 1/8 x 4 in.  | 1 1/2 x 2 3/8 in. | 50            | 8 boxes       |
|  | 7238002   | A6254 | 3 1/8 x 6 in.  | 1 1/2 x 4 3/8 in. | 50            | 8 boxes       |
|  | 7238003   | A6254 | 4 x 7 7/8 in.  | 2 x 6 3/8 in.     | 50            | 8 boxes       |
|  | 7238004   | A6254 | 4 x 10 in.     | 2 x 8 1/8 in.     | 50            | 8 boxes       |
|  | 7238005   | A6255 | 4 x 11 7/8 in. | 2 x 9 3/4 in.     | 50            | 8 boxes       |
|  | 7238006   | A6255 | 4 x 13 5/8 in. | 2 x 11 1/4 in.    | 50            | 8 boxes       |

| Leukomed® T   | Order No. | HCPCS | Dressing Size     | Wound Pad Size | Dressings/Box | Shipping Unit |
|---|-----------|-------|-------------------|----------------|---------------|---------------|
|  | 7238100   | A6257 | 3 x 2 in.         | N/A            | 50            | 10 boxes      |
|   | 7238101   | A6257 | 3 1/8 x 4 in.     | N/A            | 50            | 10 boxes      |
|   | 7238108   | A6258 | 4 x 4 3/4 in.     | N/A            | 50            | 8 boxes       |
|   | 7238102   | A6258 | 4 1/4 x 5 1/2 in. | N/A            | 50            | 10 boxes      |
|   | 7238109   | A6258 | 6 x 7 7/8 in.     | N/A            | 50            | 6 boxes       |
|   | 7238110   | A6259 | 6 x 10 in.        | N/A            | 50            | 6 boxes       |

| Leukomed® T plus  | Order No. | HCPCS | Dressing Size  | Wound Pad Size    | Dressings/Box | Shipping Unit |
|---|-----------|-------|----------------|-------------------|---------------|---------------|
|  | 7238200   | A6203 | 2 x 3 in.      | 1 1/8 x 1 1/2 in. | 50            | 8 boxes       |
|   | 7238201   | A6203 | 3 1/8 x 4 in.  | 1 1/2 x 2 3/8 in. | 50            | 8 boxes       |
|   | 7238202   | A6203 | 3 1/8 x 6 in.  | 1 1/2 x 4 3/8 in. | 50            | 8 boxes       |
|   | 7238214   | A6203 | 4 x 7 7/8 in.  | 2 x 6 3/8 in.     | 50            | 8 boxes       |
|   | 7238203   | A6204 | 4 x 10 in.     | 2 x 8 1/8 in.     | 50            | 6 boxes       |
|   | 7238204   | A6204 | 4 x 11 7/8 in. | 2 x 9 3/4 in.     | 50            | 5 boxes       |
|   | 7238205   | A6204 | 4 x 13 5/8 in. | 2 x 11 1/4 in.    | 50            | 4 boxes       |

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
 Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com  
 To order toll-free: BSN medical (+1) 800 552 1157

53110 R2

© 2013 BSN medical Inc.




REV 01/13

## Leukomed® Sterile wound protection



Well Protected –  
Free to Breathe





## Leukomed® Composite Dressing

- Highly Absorbent Dressing




## Leukomed® T Transparent Film Dressing

- Waterproof
- Bacterial Barrier



## Leukomed® T Plus Post-Op Dressing

- High Absorbent Pad
- Waterproof Film
- Bacterial Barrier



Substrate:

Soft, conformable  
non-woven



Substrate:

Stretchable, transparent  
polyurethane film

Pattern coat technology  
–ensures high moisture vapor permeability



Substrate:

Stretchable, transparent  
polyurethane film

Pattern coat technology  
–ensures high moisture vapor permeability




**Absorbs 12x its weight**  
– allows for less frequent  
dressing changes




**Non-adherent wound pad**  
– prevents sticking and  
interruption of the  
healing process



**Permeable to air**  
– reduces risk of maceration



**Waterproof and  
bacteria proof**  
– ultra-thin transparent  
polyurethane film protects  
the wound site




**Permeable to  
water vapor and air**  
– reduces risk of maceration



**Easy-to-use**  
– stabilizes the dressing during  
application without tension  
or creases



**Absorbs 12x its weight**  
– allows for less frequent  
dressing changes



**Waterproof and  
bacteria proof**  
– ultra-thin transparent  
polyurethane film protects  
the wound site



**Non-adherent wound pad**  
– prevents sticking and  
interruption of the  
healing process